

Grossing Station

High grade stainless steel 304 is used to construct the Cominus Grossing Work Station, which has backdraft ventilation and continuous air flow throughout the workspace. There are two sprinkler systems there.

Features



STAINLESS STEEL PERFORATED WORKTABLE

Stainless Steel removable perforated worktable.



HAND SPRAY ASSEMBLY

Flexible hand spray assembly helpful for user to clean the entire work area.



ILLUMINATION WITH UV LIGHT & SPOT LIGHT

Inbuilt chamber illumination for complete work area. It has an additional spot light with movable hose to illuminate a particular area.



MAGNETIC TOOL HOLDER

Magnetic tool holder to keep all the necessary tools in place.



LASER INTEGRATED RULER

Laser integrated ruler for easy measurements. Optional removable ruler is also available.




UV DISINFECTION LIGHT

UV disinfection light for disinfecting the entire work area.



Technical Specifications:-

- **Construction complete AISI 304**
- **Workbench for 2 person**
- **Dimensions 11*3.5*7 (W*D*H)**
- **Perforated Platform**
- **2 Nos Sink (Left &Right)**
- **Accessories Fitting:-Hand Towel Rod, tissue paper roll, magnetic tool holder, Steel ruler, Writing pad**
- **2 Nos Shelf**
- **Illumination- LED and UV Light**
- **Magnifying Glass**
- **Level Adjuster**
- **Ventillation With Flexible Hose**



Check for the latest cominus products

<https://cominus.in/>

Office:- D-132 Sec 4 Bawana

DSIIDC Complex Delhi-110039

Landline:- 011-47553098

Mob:- 9027822747, 9667529087 Email:- sales@cominus.in