

## Incubator Shaker

### A Reliable mixing, shaking and temperature control

With wide temperature range 15 to 80 °C and a connectable circular movement, “Cominus” shaking incubators are suitable for applications such as homogenization, incubation, fermentation or chemical and biochemical reactions.



### Key Features

- Perfect appearance with display area for temperature, speed and timing. With the enlarged digital display and clear symbols on the monitor, you can observe from a greater distance
- Sliding black window, easy to push and pull for dark culture (Optional)
- Double glass doors ensure excellent insulation and security
- Brushed full stainless steel rounded corners of the integrated cavity, beautiful and easy to clean
- Machine operation is nearly silent, multi-unit stacked high-speed operation without abnormal vibration
- Waterproof fan without heat, significantly reducing background heat and saving energy
- Flexible placement, stackable, effective in saving lab space
- Multi-safety design for operator and sample safety



## Configuration List

Incubator Shaker	1
Shaking Tray	1
Static Shelf	1
Power Cord	1
Product Manual	1
Test Report	1

## Technical Details

<b>Model No</b>	<b>CM-IS-1818</b>
Operating Temperature	15 to 80°C
Temperature Accuracy	0.5°C
Temperature Uniformity	0.3°C
Speed Range	15 to 400
Motor	BLDC
Orbit	25 mm
Timed & Cont. Operation	Upto 99 hrs
Audio Visual Alarm	Yes
Flask Fitting	
125 ml	8
250 ml	4
1000 ml	4
2000 ml	2
Temperature Sensor	PT100
Platform Size	18" x 18"
Programmable Features	Yes
Slide Out Platform	Yes
CFC Free Refrigerant	Yes
Stackable	Yes
Material of Construction	Internal AISI304 Outer Body MS Powder Coated
Door	Double Door Internal Glass External Metal

Check for the latest cominus products

<https://cominus.in/>

**Office:- D-132 Sec 4 Bawana**

**DSIIDC Complex Delhi-110039**

**Landline:- 011-47553098**

**Mob:- 9027822747, 9667529087 Email:- sales@cominus.in**