

Our bench-top orbital shakers are made to produce uniform mixing for productive cultivations. They can support any type of laboratory vessel, including flasks, tubes, and microplates. They easily fit into the majority of incubators and fume hoods thanks to their compact form.



Microprocessor PID



Continuous operation



Intuitive LCD controller displays time & RPM control.



Built-in power interruption protection function



Easy-set digital timer (1 - 9999min)



Over current protection



Accessories



Universal Tray



Lotus Clamp

Technical Specifications

Model No		<i>CM-OS-12</i>
Control System		PID Microcontroller
Display		LCD
Platform Size	Inches	12x12x12
Max Capacity		2000 ml
Maximum Load	kgs	4
Net Weight	Kgs	15
Speed Range	RPM	50-300
Amplitude	MM	10/20
Timer	Min	1-9999
Power		220 V AC 50 Hz
Warranty		1 Year

Check for the latest cominus products

<https://cominus.in/>

**Office:- D-132 Sec 4 Bawana
DSIIDC Complex Delhi-110039**

Landline:- 011-47553098

Mob:- 9027822747, 9667529087

Email:- sales@cominus.in