



Hot Air Oven

Discover our range of high-quality hot air ovens designed to meet your laboratory and industrial heating needs. Cominus is committed to providing reliable and efficient equipment, and our hot air ovens are no exception. With a capacity range from 90 to 400 Liters and a temperature range from ambient to 350 degrees Celsius, our ovens offer precise temperature control and uniform heat distribution for various applications.

Features

- The 90 to 400-liters capacity ensures sample space for accommodating large quantities of samples, making it ideal for laboratories and industrial settings with high-volume requirements.
- The oven's advanced air circulation system ensures uniform heat distribution throughout the chamber, minimizing temperature variations and ensuring consistent heating across all samples.
- Cominus Scientific hot air ovens are designed with safety in mind. They feature built-in over-temperature protection and a safety thermostat to prevent overheating and ensure the well-being of your samples and equipment.
- The user-friendly interface and intuitive controls make operating our hot air ovens a breeze. With adjustable shelves and a removable chamber for easy cleaning, maintenance is convenient and hassle-free.





Technical Specifications: -

Model	CM-HAO-90	CM-HAO-200	CM-HAO-300	CM-HAO-400
Capacity (Litres)	90	200	300	400
Temp Range (°C)	Amb +5 to 350 °C	Amb +5 to 350 °C	Amb +5 to 350 °C	Amb +5 to 350 °C
Temp Accuracy (At 37 °C)	±0.1	±0.1	±0.1	±0.1
Temp Uniformity (°C)	±0.1	±0.1	±0.1	±0.1
Glass View Door	Yes	Yes	Yes	Yes
Digital Door	Yes	Yes	Yes	Yes
Microcontroller Based	Yes	Yes	Yes	Yes
No. of Shelves	2	2	2	2
Inner Dimension	450 x 450 x 450	550 x 550 x 800	730 x 550 x 750	760 x 660 x 800
Overall Dimension	640 x 775 x 1210	740 x 875 x 1380	920 x 875 x 1410	950 x 985 x 1560
MOC internal	SS304	SS304	SS304	SS304
MOC Body	MS	MS	MS	MS
Insulation	Ceramic Wool	Ceramic Wool	Ceramic Wool	Ceramic Wool
Castor Wheel	Yes	Yes	Yes	Yes

Check for the latest cominus products

<https://cominus.in/>

Office:- D-132 Sec 4 Bawana

DSI IDC Complex Delhi-110039

Landline:- 011-47553098

Mob:- 9027822747, 9667529087

Email:- sales@cominus.in