



Vacuum Oven

We are pleased to introduce you to our high-quality vacuum ovens designed to meet your specific needs. Our vacuum ovens are engineered for various applications that require precise temperature control and a controlled environment.

Operating at a maximum temperature of 200 degrees Celsius, our vacuum ovens offer exceptional temperature uniformity and stability throughout the chamber. This ensures accurate and reliable heating for your critical processes.

With a generous capacity of 90 Liters, our vacuum ovens provide ample space for accommodating a wide range of samples, parts, or materials. Whether you need to dry, cure, or degas your products, our ovens offer sufficient room to handle larger volumes while maintaining optimal temperature conditions.

These vacuum ovens are equipped with advanced features, such as a user-friendly control panel for easy operation, digital temperature display for precise monitoring, and programmable settings for automated temperature profiles. The vacuum system ensures a controlled atmosphere inside the chamber, allowing for efficient heat transfer and elimination of unwanted contaminants.





Features

- Maximum operating temperature: 200 degrees Celsius
- Generous capacity: 90 liters
- Precise temperature control for accurate heating
- Excellent temperature uniformity and stability
- Vacuum system for a controlled atmosphere
- Digital temperature display for easy monitoring
- Programmable settings for automated temperature profiles
- User-friendly control panel for easy operation
- Durable construction for long-lasting performance
- Safety features including over-temperature protection and door interlocks
- Versatile and suitable for various applications in research, industry, and scientific settings

Technical Specifications: -

Model	CM-VO-90
Capacity (Litres)	90
Temp Range (°C)	Amb +5 to 200 °C
Temp Accuracy (At 37 °C)	±0.5
Temp Uniformity (°C)	±0.5
Glass View Door	Yes
Digital Door	Yes
Microcontroller Based	Yes
No. of Shelves	8
Inner Dimension (mm)	450 x 450 x 450
Overall Dimension (mm)	640 x 775 x 1210
MOC internal	SS304
MOC Body	SS
Vacuum Range	0 – 0.1 MPa
Castor Wheel	Yes

Check for the latest cominus products

<https://cominus.in/>

Office:- D-132 Sec 4 Bawana

DSI IDC Complex Delhi-110039

Landline:- 011-47553098

Mob:- 9027822747, 9667529087

Email:- sales@cominus.in