

Slide Cabinet

Introducing the Slide Cabinet!!

Cominus Scientific is here with Slide Cabinets for the storage of slides prepared by the painstaking process. This slide cabinet is very convenient as its slides are 75 x 25 mm and it is made as per industry standards.

Slides are made up of CRCA sheets with powder coated finish after 7-tank phosphating process which provides resistant to scratch and make it rust free even user uses it roughly.

At Cominus Scientific, we talk about 'Quality and Precision' and our target is to add it in our every product so as in the Slide Cabinets.

Our Slide Cabinet has the better ability to store slides and avoid the fungal growth which you can not see in the cabinets of wood & sun-mica and the final touch of powder coating makes it a quality and industry ready product.

Every Drawer in the slide cabinet moves smoothly in their slot and it easily removable so user can easily lift the slides from the drawers. This slide cabinet is made in such a way that it can utilize optimum space and for security we have lockable doors.

Technical Specifications: -

Model	Nos of Slides	Nos of Drawers
CM-SC-25	25000	80
CM-SC-50	50000	160
CM-SC-100	100000	336

COMINUS

SCIENTIFIC



Check for the latest cominus products

<https://cominus.in/>

Office:- D-132 Sec 4 Bawana

DSIIDC Complex Delhi-110039

Landline:- 011-47553098

Mob:- 9027822747, 9667529087

Email:- sales@cominus.in